

29th Congress of the European Academy of Allergy and Clinical Immunology London, 5th to 9th June, 2010

Differential allergen sensitization patterns in common Cypress (Cupressus sempervirens) pollen allergy: identification of a novel alkaline allergen

Shahali Y¹, Sutra JP¹, Haddad I², Vinh J², Mari A³, Peltre G¹, Charpin D⁴, Sénéchal H^{1, 5}, Poncet P^{1,6} 1- ESPCI ParisTech, PECSA, UMR 7195 CNRS, Paris, France; 2- ESPCI ParisTech, USR 3149 CNRS, Paris, France; 3- Centre for Clinical and Experimental Allergology, IDIIRCCS, Rome, Italy4- Hôpital Nord, Marseille, France ; 5- Inserm, Paris, France; 6- Pasteur Institute, Paris, France



sempervirens pollen extracts. The present study led to the revelation of at least a novel basic allergen as well as differential allergen sensitization patterns. These novel IgE reactive components may subsequently be applied to expand the panel of well-defined cypress pollen molecules for a more efficient allergen-based diagnosis and therapy.



Mediterranean basin and pollen represents the principal cause of late winter-early spring allergy. Up to now, two Cs allergens have been well characterized:

Cup s 1, a 45 kDa protein within the pectase-lyase family and **Cup s 3**, a thaumatin-like protein of **34** kDa .


Till now, most of the reported investigations on pollen molecular allergens were based on the extraction and characterization of water soluble fractions and very few attempts have been made to evaluate the allergenic potency of water insoluble pollen constituents.

Cypress pollens are particularly liable to desiccation during their transport, facilitating the exudation and accessibility of both hydrophobic and hydrophilic intra-pollinic materials to the environment. In order to bring to light the potential role of most insoluble fractions in the cypress pollen allergy, we choose a proteomic approach based on the extraction and solubilisation of both hydrophilic and hydrophobic proteins by using some chaotropes and detergents.

Kinetic of allergen liberation in cypress pollen grains Dry conditions Wet conditions

In the usual case, the passage through the pollen exine encounters a high and virtually constant resistance. However in cypress pollen, the fragility of the exine facilitates the exudation and accessibility of hydrophobic intrapollinic materials to the environment.

Efficiency of the protein extraction



1-D (SDS-PAGE) immunoblotting using SDS-extracts: Differential IgE binding patterns in *C. sempervirens* pollen allergy





Figure1. IgE reactivity to C. sempervirens SDS-pollen extract. Mr represents the marker; Cs, C. sempervirens pollen extract; lanes 1-24, screening of 24 patient's sera by Western blotting for specific IgE antibodies to C. sempervirens SDS- pollen extracts and Co, the negative control.






Figure 2. Silver staining of 2-DE separation (G) and typical 2-DE IgEbinding spectra of C. sempervirens pollen proteins (A, B, C). Molecular masses are expressed in kilo-Daltons (kDa); pl corresponds to isoelectric points. Representative IgE binding patterns A, B and C correspond respectively to 2-DE immunoblotting allergen profiles A, B and C shown in the Fig.1.





Figure 3. Comparative 2D-PAGE using H₂O (a, c) and PBS (b, d) extracts. Immunoprints (e, f) were performed with a same intermediate serum reacting to both 14 kDa allergen and high molecular weight proteins. This picture clearly demonstrates that different extraction methods result to different 2-DE allergen profiles. Aqueous extraction was very selective and essentially gave access to several isoforms of a 43 kDa protein (frame 1), while the PBS allowed the extraction of some additional basic allergenic spots of 43, 60 kDa (frame 2) as well as three basic isoforms of about 14 kDa.



MS/MS results						
Protein Name	Accession No.	Рер				

Frame	Spot (s)	Protein Name	Accession No.	Peptide count	Mascot score		
1	1	Cup s 1 pollen allergen precursor (pectate lyase family)	gi 8101715	17	551		
2	2	*Putative allergen Cup a 2 variant 1 [Hesperocyparis]- Polygalacturonase	gi 11819795	8	120		
2	3	*Polygalacturonase OS= <i>Juniperus ashei -</i> GN=JNA2	PGLR2_JUNAS	8	133		
3	4	Histone H4	B6TJ90_MAIZE	6	308,21		
3	5	Histone H4	B6TJ90_MAIZE	8	437,43		
3	6	🗡 Histone H4	B6TJ90_MAIZE	4	210,31		
Candidate, has yet to be confirmed							

allergen presents in Cupressaceae pollen grains isolated and reported for the first time in C. sempervirens